

GLAM® Ceramic Brackets

Advantages at a glance

- Very discreet – matt, translucent ceramic
- Permanent aesthetics – no discolouration
- Solid – made of deformation-free material
- Reliable – the forces of the archwires are conveyed optimally to the tooth
- Comfortable – rounded edges provide high wearing comfort
- Caring approach – they can be removed without effort at the end of treatment
- No allergies – made of medical ceramics




Reliability & precision made in Germany

FORESTADENT® is a company specialising in manufacturing orthodontic materials and is among the world leaders in this segment. The company was founded in 1907 as a manufacturer of jewellery and produced watches until well into the 1970's with their own movement.

Even today, FORESTADENT® benefits from its know-how in manufacturing minute high-precision parts in combination with advanced surface finishing. The company has been family-run for four generations. For over 100 years, FORESTADENT® has been manufacturing in its main facility in Pforzheim and meanwhile markets its products in over 80 countries.

Ask your orthodontist for advice.



GLAM® Ceramic Brackets

The elegant approach
to a beautiful smile



The elegant approach to a beautiful smile

Wearing braces is absolutely normal these days. After all, a beautiful smile becomes truly attractive with straight teeth. Straight teeth are also easier to care for and remain healthier. Wearing braces is no longer a prerogative of young people, more and more adults have started wearing them. Even celebrities occasionally appear in public wearing braces. And although it may be standard practice these days to have one's teeth corrected, not everybody wants to be seen wearing braces.

GLAM® solves the problem. The most prominent part of fixed braces are the so-called brackets which are bonded to the teeth. GLAM® Brackets are made of special medical ceramic material allowing the natural colour of the tooth to shine through. This makes your braces virtually invisible and already gives you a confident smile during treatment.

Optimal aesthetics

The material for the GLAM® Brackets is not coloured solid white, but is translucent. This means, the natural colour of your teeth can shine through giving an optimal camouflage effect. The surface of GLAM® Brackets has a matt glow instead of being glossy, which prevents any irritating light reflections. And best of all: ceramic does not discolour, it retains its aesthetic properties throughout the entire treatment period.

Comfortable and safe

The mode of action of ceramic brackets is no different than that of metal brackets. They do not lose their shape and convey the clinical forces necessary for positioning teeth optimally. The edges of the GLAM® Brackets are rounded to provide as little irritation as possible. GLAM® Brackets remain secure and fixed during the entire treatment period and can be removed with little effort at the end. Another benefit: GLAM® Ceramic Brackets do not require any metal parts which makes them suitable for patients suffering from allergies.

Perfect combination

Fixed braces do not only consist of brackets, but also include elastic archwires. These are inserted into the brackets and provide the actual forces necessary to move the teeth. When deciding in favour of high-quality brackets, one should also consider archwires of corresponding quality. The FORESTADENT® BioCosmetic® archwires are especially suited for GLAM®. They consist of an ultra-elastic metal core covered with a tooth-coloured layer of Teflon. This ensures that they remain effective and permanently white. Further information on the FORESTADENT® BioCosmetic® and other archwires is given in our separate "BioArchwires" brochure.

