

# ORTHOeasy® Pal


The new design of the thread peak offers safe and easy passage through the gingiva and corticalis.

Shark like cuts make a new pitch at every single winding step.

Transgingival conic screw neck for sealing the gingiva against infections.

Compression step and depth stop at the corticalis for higher stability through mechanical compaction of the bone.

Easy anchorage of appliances in the screw head.



## OrthoEasy Pal – the new palatal pin

The OrthoEasy System is now complemented by the OrthoEasy Pal. Pal stands for palatal, as the new pin was developed specifically for the anchorage of orthodontic appliances in the palate. For example, the Frog appliance can now be attached quicker and easier in the mouth than before. Until now, the lab abutments with the welded appliance had to be attached elaborately with a wire ligature on the pins placed in the palate. The head of the new OrthoEasy Pal is equipped with a practical inner thread so that the abutments can simply be attached with a retaining screw. A lab analog pin and an impression cap help to fabricate the appliance in the orthodontic laboratory.

## Registration

Please fill in your details below and fax to 01908 320029 or email to [info@forestadent.co.uk](mailto:info@forestadent.co.uk) to express your interest and we will call you back to provisionally hold a place. A place will only be guaranteed upon full payment\*. Alternatively call 01908 227 851 for more details.

Yes, I would like to participate in the course

13<sup>th</sup> - 14<sup>th</sup> October, 2017 - London

Name

Address

Phone

Fax

E-mail

GDC no.

Date

Signature

Practice stamp

\*Cancellation policy - All bookings cancelled at least 6 weeks before the date of the course will be entitled to a 50% refund or a full transfer to an alternative date within the next year. Any cancellations made less than 6 weeks before the course will be entitled to a full transfer to an alternative date within the next year but will not receive any refund.

# ORTHOeasy®

## Cortical Anchorage Course:

# Miniscrews & Anchorage in Modern Orthodontics



The Complete System for Cortical Anchorage

Friday 13<sup>th</sup> & Saturday 14<sup>th</sup>  
October 2017, London  
£649 + VAT



## Aims of the course

*Anchorage: The enemy or friend that we have to deal with every day during our clinical practice.*

*To recognise the need for anchorage and to know how to obtain and maintain it during treatment are the key points that lead to orthodontic success. For years, orthodontists have applied complex systems and various appliances for anchorage, sometimes only succeeding in part.*

*The absence of multiple teeth or bone loss can sometimes make orthodontic treatment extremely complex, if not impossible. The introduction of Temporary Anchorage Devices (TADs) has revolutionized the orthodontic world, simplifying otherwise complex treatments and making impossible treatments possible.*

*Like other innovative procedures, the application of TADs can be initially rejected by the clinician, who may consider it complex and unnecessary. The course is intended to illustrate the advantages and disadvantages of TADs, their features and benefits, simple clinical procedures, as well as the possible risks and complications. Our goal will be to increase the knowledge and confidence of the clinical application of TADs, through the analysis of different clinical situations, as well as by applying TADs on blocks of synthetic bone.*



Colour coded anodising of the surface.

## Speaker



### Giorgio Iodice

*Giorgio Iodice obtained his degree in Dentistry, his specialist degree in Orthodontics and his PhD in Oral Sciences from the University of Naples Federico II. Giorgio received the Certificate of Excellence in Orthodontics in Italy (Italian Board of Orthodontics) and the E.B.O. (European*

*Board of Orthodontics). Since 2013, Giorgio has been a member of the IBO examining committee and of the SIDO Model Display examining committee since 2014. He is assistant Clinical Professor at the University of Naples Federico II and an active member of SIDO, EOS, SIBOS, AIDOr and SIDA.*

*His main clinical and research interests are focused on the characteristics of treating adult patients, temporary anchorage devices (TADs) and the relationship between occlusion and temporomandibular disorders.*

*Giorgio has written several publications in national and international journals as well as a valued speaker at many congresses and universities.*



*The Frog Appliance is suitable for effective molar distalisation without tipping and can be anchored palatally with the aid of OrthoEasy® Pal. An advanced version of this commonly used appliance is now available under the name Frog II.*

## Course program

### 1<sup>st</sup> day

#### Theory:

- Anchorage in orthodontics
- Skeletal anchorage in orthodontics
- Clinical indications for use of TADs in orthodontics
- TADs Features and benefits
- Choosing the ideal site and the type of TADs
- Clinical Procedures
- Risks and Complications
- Clinical show cases as examples of the different clinical situations
- Auxiliary systems on TADs

#### Practice:

- Clinical show cases as examples of the different clinical situations
- TADS insertion on blocks of synthetic bone.

### 2<sup>nd</sup> day

- Systems with skeletal palatal anchorage: Ortho Easy Pal
- Clinical show cases as examples of the different clinical situations



OrthoEasy® Pal – the new palatal pin.