




Track[®] Thermoforming foils

Faster is better



www.forestadent.com

FORESTADENT[®]
GERMAN PRECISION IN ORTHODONTICS

TRACK® Thermoforming foils



Fast processing due to integrated separating foil!

Track® A and B cut down on several processing steps: the use of separating liquid or a separating foil and manual preparation is not required because Track® A and B have an integrated separating foil.



No bond enhancing agent required!

Track® A and Track® B thermoforming foils do not require bond enhancing medium to bond to acrylic (Forestacryl).



Superior thermobonding characteristics!

Better thermobonding when separately processed parts of an appliance are joint together i.e. snoring appliances, positioners and mouthguards. The improved performance during thermobonding reflects the use of superior materials in the manufacturing process of Track® E.

Ø 125 mm, transparent

Order number:	thickness	contents		description	Application	Bonding to acrylic
408-0105	0,5 mm	100		Track® A Hard (PETG)	invisible retainers, standard retainers, correction splints, applications like Essix A.	yes
408-0106	0,6 mm	100				
408-0108	0,8 mm	100				
408-0110	1,0 mm	100				
408-0115	1,5 mm	50				
408-0120	2,0 mm	50				
408-0220	2,0 mm	30		Track® B Double layer (TPU/PET-G)	comfortable retainers with optional correction (if used with set-up)	yes (on hard side)
408-0305	0,5 mm	100		Track® C Soft (PE)	transfer trays, etching trays, applications like Essix C	no
408-0310	1,0 mm	100				
408-0515	1,5 mm	50		Track® E Soft, rubbery (EVA)	positioners, bracket transfer trays, bracket guards, positioners, can be used to duplicate study models	no
408-0520	2,0 mm	50				
408-0530	3,0 mm	30				
408-0610	1,0 mm	100		Track® bleach (EVA)	bleaching trays, applicating tray for flouride and medication	no

TRACK[®]-V vacuum thermoforming unit












- ① Medium-wave infrared heater with cover
- ② Foil reception handle with cavity for good grip
- ③ Form pot for reception of model, granules pot and model plate (Occluform-3, 188 580)
- ④ Main switch
- ⑤ Granules collecting device

- ⑥ Foil reception with foil securing ring
- ⑦ Mains connection (back)
- ⑧ Sensor window for temperature registration
- ⑨ Column for movement of oil reception
- ⑩ Touch panel to operate all functions

Order number:		contents	description
900-2000		1	Track[®]-V vacuum thermoforming unit, 230 V
900-2110		1	Track[®]-V vacuum thermoforming unit, 100 V, Japan
900-2115		1	Track[®]-V vacuum thermoforming unit, 115 V, USA

TRACK®-V auxiliaries

Order number:	contents		description
408-0705	Cover template , ø 0,5 x 100 mm (8 pcs.)		1
408-0706	Occluform for Track®-V		1
900-2001	Filling granules , refill package 1,3 kg		1
900-2002	Blocking out plaster , 400 g		1
900-2003	Blocking out wax , transparent, 50 g		1
900-2012	Brush for granules , incl. magnet		1
900-2016	Sealing for foil reception		1
900-2017	Sealing for form pot		1
900-2018	Finishing set Track® for cutting and handling tracks		1