

Currently validated printers

DLP technology

Resin type	SprintRay Pro 95	SprintRay Pro 55	Carima IM2	Carima IMD (IM3)	Asiga Pro 4 k	Asiga MAX UV
TC-85 DAW & TC-85 DAC Direct Aligners	✓	✓	✓	✓	✓	✓
TC-80 DP Permanent C&B	✓	✓	✓	✓	✓	✓
SG-100 Surgical Guide	✓	✓	✓	✓	✓	✓
THD Denture Base	✓	✓	✓	✓	✓	✓
TFDH Flexible Denture	✓	✓	✓	✓	✓	✓
S-100M Model	✓	✓	✓	✓	✓	✓
SC-130 Castable	✓	✓	✓	✓	✓	✓
TE-Series Mouth & Night Guard	✓	✓	✓	✓	✓	✓

Currently validated printers

LCD technology

Resin type	Uniz SLASH 2	Uniz NBEE	Shinwon Dental Kave	Ackuretta Sol	Cybermed OnDemand 3D	Phrozen Shuffle	Sindoh SD	Ray Ram600
TC-85 DAW & TC-85 DAC Direct Aligners	✓	✓	✓	✓	✓	✓	✓	✓
TC-80 DP Permanent C&B	✓	✓	✓	✓	✓	✓	✓	✓
SG-100 Surgical Guide	✓	✓	✓	✓	✓	✓	✓	✓
THD Denture Base	✓	✓	✓	✓	✓	✓	✓	✓
TFDH Flexible Denture	✓	✓	✓	✓	✓	✓	✓	✓
S-100M Model	✓	✓	✓	✓	✓	✓	✓	✓
SC-130 Castable	✓	✓	✓	✓	✓	✓	✓	✓
TE-Series Mouth & Night Guard	✓	✓	✓	✓	✓	✓	✓	✓

Currently validated printers

SLA technology

Resin type	Sindoh A1+
TC-85 DAW & TC-85 DAC Direct Aligners	✓
TC-80 DP Permanent C&B	✓
SG-100 Surgical Guide	✓
THD Denture Base	✓
TFDH Flexible Denture	✓
S-100M Model	✓
SC-130 Castable	✓
TE-Series Mouth & Night Guard	✓