




# *Tulip buccal tubes*



*Faster ligating – more patient comfort*



# Tulip buccal tubes

The Tulip buccal tubes owe their name to the tulip-shaped entrance funnel. It has been dimensioned especially large to allow quick and easy insertion of the archwires. Your patients will be delighted by the extremely flat design which guarantees a high level of comfort.

Tulip buccal tubes are available as single buccal tubes with hooks on Large Pad or on the extra big Big Foot Pad.



Maxillary	Tooth	Torque	Angulation	distal Offset	Slot	on Large Pad		on Big Foot Pad	
						right	left	right	left
76   67	0°	-	-	0°	.018" x .028"	743L0714	742L0714	743B0714	742B0714
					.022" x .028"	743L0712	742L0712	743B0712	742B0712
76   67	-10°	-	-	0°	.018" x .028"	743L0724	742L0724	743B0724	742B0724
					.022" x .028"	743L0722	742L0722	743B0722	742B0722
76   67	-10°	-	-	6°	.018" x .028"	743L0734	742L0734	743B0734	742B0734
					.022" x .028"	743L0732	742L0732	743B0732	742B0732
76   67	-14°	-	-	10°	.018" x .028"	743L0744	742L0744	743B0744	742B0744
					.022" x .028"	743L0742	742L0742	743B0742	742B0742
76   67	-14°	-	-	14°	.018" x .028"	743L0754	742L0754	743B0754	742B0754
					.022" x .028"	743L0752	742L0752	743B0752	742B0752

Mandibular	Tooth	Torque	Angulation	distal Offset	Slot	on Large Pad		on Big Foot Pad	
						right	left	right	left
7   7	-10°	-	-	0°	.018" x .028"	742T0724	743T0724	-	-
					.022" x .028"	742T0722	743T0722	-	-
76   67	-20°	-	-	0°	.018" x .028"	743L1714	742L1714	743B1714	742B1714
					.022" x .028"	743L1712	742L1712	743B1712	742B1712
76   67	-25°	-	-	6°	.018" x .028"	743L1724	742L1724	743B1724	742B1724
					.022" x .028"	743L1722	742L1722	743B1722	742B1722
76   67	-30°	-	-	6°	.018" x .028"	743L1734	742L1734	743B1734	742B1734
					.022" x .028"	743L1732	742L1732	743B1732	742B1732

## FORESTADENT Ltd.

Unit 1 · Crossinglands Business Park · Salford Road  
Aspley Guise · MK17 8HZ · Buckinghamshire · Great Britain  
Phone: +44 1908 227 851 · info@forestadent.co.uk

## FORESTADENT Bernhard Förster GmbH

Westliche Karl-Friedrich-Straße 151 · D-75172 Pforzheim (Germany)  
Phone: +49 (0) 7231 459-0 · Fax: +49 (0) 7231 459-102  
info@forestadent.com · www.forestadent.com