



New!

Tulip buccal tubes

Faster ligating – more patient comfort



www.forestadent.com

FORESTADENT®
GERMAN PRECISION IN ORTHODONTICS

Tulip buccal tubes

New!

Tulip buccal tubes

The new Tulip buccal tubes owe their name to the tulip-shaped entrance funnel. It has been dimensioned especially large to allow quick and easy insertion of the archwires. Your patients will be delighted by the extremely flat design which guarantees a high level of comfort.



Maxillary

Tooth	Torque	Angulation	distal Offset	Slot	on Large Pad		on Big Foot Pad	
					right	left	right	left
$\frac{76}{67}$	0°	–	0°	.018" x .025"	743L0714	742L0714	743B0714	742B0714
$\frac{76}{67}$	0°	–	0°	.022" x .025"	743L0712	742L0712	743B0712	742B0712
$\frac{76}{67}$	-10°	–	0°	.018" x .025"	743L0724	742L0724	743B0724	742B0724
$\frac{76}{67}$	-10°	–	0°	.022" x .025"	743L0722	742L0722	743B0722	742B0722
$\frac{76}{67}$	-10°	–	6°	.018" x .025"	743L0734	742L0734	743B0734	742B0734
$\frac{76}{67}$	-10°	–	6°	.022" x .025"	743L0732	742L0732	743B0732	742B0732
$\frac{76}{67}$	-14°	–	10°	.018" x .025"	743L0744	742L0744	743B0744	742B0744
$\frac{76}{67}$	-14°	–	10°	.022" x .025"	743L0742	742L0742	743B0742	742B0742
$\frac{76}{67}$	-14°	–	14°	.018" x .025"	743L0754	742L0754	743B0754	742B0754
$\frac{76}{67}$	-14°	–	14°	.022" x .025"	743L0752	742L0752	743B0752	742B0752

Mandibular

Tooth	Torque	Angulation	distal Offset	Slot	on Large Pad		on Big Foot Pad	
					right	left	right	left
$\frac{76}{67}$	-20°	–	0°	.018" x .025"	743L1714	742L1714	743B1714	742B1714
$\frac{76}{67}$	-20°	–	0°	.022" x .025"	743L1712	742L1712	743B1712	742B1712
$\frac{76}{67}$	-25°	–	6°	.018" x .025"	743L1724	742L1724	743B1724	742B1724
$\frac{76}{67}$	-25°	–	6°	.022" x .025"	743L1722	742L1722	743B1722	742B1722
$\frac{76}{67}$	-30°	–	6°	.018" x .025"	743L1734	742L1734	743B1734	742B1734
$\frac{76}{67}$	-30°	–	6°	.022" x .025"	743L1732	742L1732	743B1732	742B1732