

IT'S ALL ABOUT THE CLIP...

BioQuick®


The 5th Generation

www.forestadent.com

FORESTADENT®
GERMAN PRECISION IN ORTHODONTICS

BioQuick®

It's all about the Clip...




The current Clip has still the original width but has been optimized and extended in the heavily stressed area. This allows it to absorb more force and is less prone to deformation.


The advantages:

- The larger Clip opening allows for easier opening with tools such as explorers or scalars. In addition, we have developed a new instrument for even safer opening.
- The application point for the force required for opening has been optimized. This means that less force is needed and at the same time there is less strain on the material.
- The new Clip design allows a larger usable slot width and hence better rotation and angulation control, while at the same time the bracket has become smaller.
- There are four Clip variants for different bracket angulations: the bracket is now opened simply and intuitively along the longitudinal axis.

The new clip is opened in the angulation direction:



Midline marking remains in open state correct:



All ordering information on BioQuick® in the respective Prescription.



IT'S ALL ABOUT THE CLIP...